

PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT FOR FUTURE USE IN CASE YOU HAVE A PROBLEM.

BAUSCH & LOMB

Sensitive Eyes® Plus Saline Solution

CONTENTS: Sensitive Eyes Plus Saline Solution is a sterile, isotonic, buffered solution that contains boric acid, sodium borate, potassium chloride, sodium chloride, preserved with polyaminopropyl biguanide (0.0003%) and edetate disodium (0.025%).

ACTIONS: Sensitive Eyes Plus Saline Solution removes loosened debris and traces of daily cleaner when used as a rinse after cleaning. Also removes existing calcium deposits and prevents new deposits from accumulating when used with heat disinfection. The saline solution may also be used in the lens carrying case during heat (thermal) disinfection. It can be used as a diluent for enzyme tablets, to store lenses, and as a final (pre-inserting) rinse after chemical (not heat) and hydrogen peroxide disinfection.

INDICATIONS (Uses): Sensitive Eyes Plus Saline Solution is a sterile, preserved solution for rinsing after cleaning and before lens insertion with heat (thermal), chemical (not heat) and hydrogen peroxide systems; heat disinfection; dilution of enzyme tablets and storage of soft contact lenses as recommended by your eye care practitioner.

CONTRAINDICATIONS (Reasons not to use): If you are allergic to any ingredient in Sensitive Eyes Plus Saline Solution, do not use this product.

PRECAUTIONS:

- Never re-use this solution.
- Keep the bottle tightly closed when not in use.
- Store at room temperature.
- Use before the expiration date marked on the carton and bottle.
- Do not alternate between heat (thermal) and chemical (not heat) methods of disinfection.
- After re-inserting your lenses, always empty your storage case, rinse with Sensitive Eyes Plus Saline Solution, and allow to air dry.

WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE. It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision.

Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight.

Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care practitioner. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement.

Studies have also shown that smokers have a higher incidence of adverse reactions.

If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner.

It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

To avoid contamination, do not touch tip of bottle to any surface. Replace cap after using.

ADVERSE REACTIONS (Problems and what to do): The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

- **IMMEDIATELY REMOVE YOUR LENSES.**
- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, **DO NOT** put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.
- If the problem continues, **IMMEDIATELY** remove the lens and consult your eye care practitioner.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

GOOD LENS CARE PRACTICES:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse lens case with fresh rinsing solution and allow to air dry.

DIRECTIONS:

Rinsing After Cleaning

- Clean your lenses with a daily cleaner such as **Bausch & Lomb® Sensitive Eyes Daily Cleaner** or with **ReNu® Multi-Purpose Solution** as instructed by your eye care practitioner.
- After cleaning, rinse lenses thoroughly with **Sensitive Eyes Plus Saline Solution** by holding the lens between the forefinger and thumb or placing it in the palm of the hand. Direct a steady stream of saline onto the lens for 10 seconds.

Heat (thermal) Lens Disinfection

- Fill the chambers of your lens carrying case 2/3 full with fresh **Sensitive Eyes Plus Saline Solution**.
- Place the cleaned lens into the correct chamber.
- Tightly close the cap of the chamber.
- Place the lens carrying case into the disinfecting unit and follow the operating instructions for the unit.
- After disinfection, you may store your lenses in the unopened case until ready to put on your eyes.

Emergency (Alternate) Method For Heat (thermal) Disinfection

If your heat disinfection unit is not available, place the tightly closed lens carrying case which contains the lenses into a pan of already boiling water. Leave the closed lens case in the pan of boiling water for at least ten minutes. (Above an altitude of 7,000 feet, boil for at least 15 minutes.) Be careful not to allow the water to boil away. Remove the pan from the heat and allow it to cool for 30 minutes to complete the disinfection of the lenses.

NOTE: Use of heat disinfection unit should be resumed as soon as possible.

Chemical Disinfection (Including hydrogen peroxide systems)

- Clean, rinse and disinfect your lenses as directed by your eye care practitioner.
- Before inserting your lenses, rinse them thoroughly with **Sensitive Eyes Plus Saline Solution**.
- Lenses that are chemically disinfected may absorb ingredients from the disinfecting solution which may cause irritation to the eye. A 20 second rinse may reduce the potential for irritation.

Enzyme (Weekly) Cleaning

- Clean and rinse lenses as directed by your eye care practitioner.
- **Sensitive Eyes Plus Saline Solution** may be used to dilute enzymatic cleaning tablets such as **Sensitive Eyes Enzymatic Cleaner**, **Bausch & Lomb ReNu Thermal** or **Effervescent Enzymatic Contact Lens Cleaner**.

HOW SUPPLIED: Bausch & Lomb Sensitive Eyes Plus Saline Solution is available for purchase in 4 Fl. Oz. (118ml) and 12 Fl. Oz. (355ml) plastic bottles. Bottles and cartons are marked with a lot number and expiration date.

PRODUCT SATISFACTION GUARANTEE
*If you are not completely satisfied with this product,
call the toll-free number on this insert.*

Toll Free Product Information
Call: 1-800-553-5340

Manufactured by:
Bausch & Lomb, Rochester, NY 14692-0450. Made in USA.