

PLEASE READ CAREFULLY AND KEEP THIS PACKAGE INSERT FOR FUTURE REFERENCE.
It is recommended that all eye care professionals review this information with their patients.

BAUSCH + LOMB

Sensitive Eyes® Saline Solution

DIRECTIONS:

Rinsing After Cleaning

- After cleaning, rinse lenses thoroughly with Sensitive Eyes® Saline Solution by holding the lens between the forefinger and thumb or placing it in the palm of the hand. Direct a steady stream of Sensitive Eyes® Saline Solution onto each side of the lens for 5 seconds, or until all traces of the cleaner are removed.

Rinsing After Chemical Disinfection (Including hydrogen peroxide systems)

- Clean, rinse and disinfect your lenses as directed by your eye care practitioner.
- Before inserting your lenses, rinse them thoroughly with Sensitive Eyes® Saline Solution.
- Lenses that are chemically disinfected may absorb ingredients from the disinfecting solution which may cause irritation to the eye. A 20 second rinse may reduce the potential for irritation.

CONTENTS:

Sensitive Eyes® Saline Solution is a sterile, isotonic, buffered solution that contains boric acid, sodium borate, potassium chloride, sodium chloride, preserved with polyaminopropyl biguanide (0.00003%) and edetate disodium (0.025%).

ACTIONS:

Bausch + Lomb Sensitive Eyes® Saline Solution removes loosened debris and traces of daily cleaner when used as a rinse after cleaning. It can also be used to rinse lens cases and as a final (pre-inserting) lens rinse after chemical (not heat) and hydrogen peroxide disinfection.

INDICATIONS (Uses):

Bausch + Lomb Sensitive Eyes® Saline Solution is a sterile, preserved solution for rinsing soft (hydrophilic) and gas permeable contact lenses after cleaning and before lens insertion with chemical (not heat) and hydrogen peroxide systems.

CONTRAINDICATIONS (Reasons not to use):

If you are allergic to any ingredient in Sensitive Eyes® Saline Solution, do not use this product.

WARNINGS: PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN CORNEAL INFECTIONS AND/OR ULCER AND LEAD TO LOSS OF VISION. It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case.

Eye problems, including corneal ulcers, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care professional. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement. Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, redness of the eye, immediately remove your lenses and promptly contact your eye care professional. It is recommended that contact lens wearers see their eye care professional twice each year or, if directed, more frequently.

To avoid contamination, do not touch tip of bottle to any surface. Replace cap after using.

WARNING: Failure to discard solution from lens case after each use or use of water to care for your lenses may lead to contamination resulting in eye injury and potential loss of vision. See accompanying instructions for additional important safety information.

Rinse Time:

- Rinse each side of the lens for 5 seconds with Sensitive Eyes® Saline Solution.

WARNING:

- You should follow the complete recommended lens rubbing and rinsing times in the product labeling to adequately clean your lenses and reduce the risk of contact lens contamination. Reduced rubbing or rinsing time may not adequately clean your lenses.
- Never use water, saline solution, or rewetting drops to disinfect your lenses. These solutions will not disinfect your lenses. Not using the recommended disinfectant can lead to severe infection, vision loss or blindness.

Lens Case Care:

- Before First Use- Rinse lens case and caps thoroughly with fresh contact lens disinfecting solution.

- Clean contact lens cases with digital rubbing with fresh, sterile disinfecting solutions/contact lens cleaner. **Never use water.** Cleaning should be followed by rinsing with fresh, sterile disinfecting solutions (**never use water**) and wiping the lens cases with fresh, clean tissue is recommended. Air-drying or recapping the lens case lids after use without any additional cleaning methods should be avoided. If air drying, be sure that no residual solution remains in the case before allowing it to air dry.
- Replace the lens case at least once every 3 months or as directed by your eye care professional.
- Contact lens cases can be a source of microbial growth.

WARNING:

- Do not store your lenses or rinse your lens case with water or any nonsterile solution. Only use fresh Sensitive Eyes Saline Solution so you do not contaminate your lenses or lens case. Use of non-sterile solution can lead to severe infection, vision loss or blindness.

Water Activity:

- Do not expose your contact lenses to water while you are wearing them.
- Ask your eye care professional for recommendations about wearing your lenses during any activity involving water.

Discard Date:

- Discard any remaining solution 90 days after opening.

WARNING:

- Using your saline solution beyond the discard date could result in contamination of the solution and can lead to severe infection, vision loss or blindness.

IMPORTANT SAFETY INFORMATION:

- Always follow the product directions for use. Failure to follow product directions may lead to vision loss.
- Visit your eye care practitioner regularly.
- Always wash and dry hands before handling lenses.
- Do not use tap water, bottled water or saliva with lenses or lens case.
- Always discard any remaining saline solution in your lens case after each use.
- Saline or rewetting drops will not disinfect your lenses.
- Always replace your solutions, lenses and lens case as directed.
- To avoid contamination, do not touch tip of container to any surface. Replace cap after using.

PRECAUTIONS:

- Never re-use this solution.
- Keep the bottle tightly closed when not in use.
- Store at room temperature.
- Use before the expiration date marked on the carton and bottle.
- Discard any remaining saline solution after 90 days.
- Do not alternate between heat (thermal) and chemical (not heat) methods of disinfection.

ADVERSE REACTIONS (Problems and what to do):

The following problems may occur: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above:

- **IMMEDIATELY REMOVE YOUR LENSES.**
- If the discomfort or problem stops, then look closely at the lens.
- If the lens is in any way damaged, DO NOT put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, or the problem stops and the lens appears undamaged, thoroughly clean, rinse and disinfect the lens, then reinsert it.

If the problem continues, IMMEDIATELY remove the lens and consult your eye care practitioner.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

When you visit your eye care professional for the assessment of a potential eye infection take your lenses, solutions and lens case to your appointment for possible culturing in order to better establish the identity of any microorganisms.

GOOD LENS CARE HYGIENE:

- Always wash and rinse your hands before you handle your lenses.
- Clean, rinse and disinfect your lenses each time you remove them.
- Always handle the same lens, the right or the left, first in order to avoid mix-ups.
- After use, always empty and rinse lens case with fresh Sensitive Eyes® Saline Solution and allow to air dry.

HOW SUPPLIED:

Bausch + Lomb Sensitive Eyes® Saline Solution is available for purchase in 12 FL. OZ. (355 mL) plastic bottles. Bottles and cartons are marked with a lot number and expiration date.

PRODUCT SATISFACTION GUARANTEE

If you are not completely satisfied with this product, call the toll-free number on this insert.

Toll Free Product Information

Call: 1-800-553-5340

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.

© Bausch & Lomb Incorporated

Manufactured by:

Bausch & Lomb Incorporated

Rochester, NY 14609